

AirGuard

Reduce medical air and oxygen misconnections

AirGuard is an accessory for medical air flowmeters that reduces gas misconnection by requiring a conscious decision with two-handed action to attach tubing to the medical air outlet.

The AirGuard provides:

- bold labelling of the air outlet
- a physical barrier that is swung aside before tubing can be connected to air

Oxygen is left unchanged and is freely accessible for emergencies.



Install AirGuards in a 400 bed hospital for the cost of 1 patient day in ICU



AirGuard is simple to operate

1 x AirGuard AUD\$25
 (flowmeter not included)
 Discounts for quantities of 100+ on application



Fitting The AirGuard is best fitted by clinical or biomedical engineering.

Quantity A 400 bed hospital with 300 air flowmeters thus needs 300 AirGuards.

About 10% more AirGuards may be needed in the 12 months after installation as more medical air flowmeters are discovered in storage.

www.fbe.com.au

The AirGuard disk is selected to fit the flowmeter.
 Flowmeters that have been fitted include:



Compact



Chrome Body

Precision Medical



Anaequip



CIG



Comweld EZI-FLOW



Therapy Equipment Diamond range