

# FBE P/L INDUROMETERS

## Viscoelastic Indurometer



### ADVANTAGES:

- Displacement with time
- Hardness of tissue
- Viscoelasticity with time
- Fluid movement
- Useable at different angles and on uneven surfaces
- Undergoing clinical trial

## Electronic Tissue Indurometer



### ADVANTAGES:

- Easy to read display
- Operate at different angles and on uneven surfaces

### DISADVANTAGES:

- Displacement only
- Snap shot in time
- No viscoelastic measurement
- Available 2011

## Mechanical Tonometer



### DISADVANTAGES:

- Displacement only
- Hardness of tissue
- Snap shot in time
- Unstable
- Hard to read inner dial
- Flat surface and vertical use only
- No longer in production